

1, 2, 3, 4
Compaction Autografting / Condensation

Maintaining Bone Bulk Results In Higher BIC

5, 6, 7, 15, 16, 17
Enhance Bone Density

Accelerates Bone Healing

8, 9, 10
Increase Residual Strain


Enhances Osteogenic Activity Through Mechanobiology

11, 12, 13, 17
Increase Implant Stability

Higher Insertion Torque & ISQ Reduces Micromotion


Versatility, It's in Our Name

Densah® Lift Protocol I




Courtesy of Dr. Salah Huwais

Densah® Lift Protocol II



Courtesy of Dr. Salah Huwais

Ridge Expansion Protocol



Courtesy of Dr. Samvel Bleyan


Ridge Expansion with Modified Ridge Split Protocol



Courtesy of Dr. Salah Huwais


Visit versah.com/versatility/ for full protocols

Immediate Implant Placement Protocol




Courtesy of Dr. Salah Huwais

Molar Septum Expansion Protocol



Courtesy of Dr. Samvel Bleyan

Guided Expansion Graft: 2-Stage Augmentation Protocol



Courtesy of Dr. Salah Huwais

Guided Surgery Protocol

Utilize our keyless guided surgery system to facilitate oral surgery.



Courtesy of Dr. Richard Martin



OSSEODENSIFICATION:
A Paradigm Shift In Implant Osteotomy Preparation
Utilizing Densah® Bur Technology



01. Todisco, M. and P. Trisi, Bone mineral density and bone histomorphometry are statistically related. *Int J Oral Maxillofac Implants*, 2005, 20(6): p. 898-904.
02. Frost HM. A brief review for orthopedic surgeons: fatigue damage (microdamage) in bone (its determinants and clinical implications). *J Orthop Sci*, 1998;3(5):272-281.
03. Kold S, et al. Bone compaction enhances fixation of hydroxyapatite-coated implants in a canine gap model. *J Biomed Mater Res B Appl Biomater*, 2005;75(1):49-55.
04. Schlegel KA, et al. Bone conditioning to enhance implant osseointegration: an experimental study in pigs. *Int J Oral Maxillofac Implants*, 2003;18(4):505-511.
05. Nkenke E, et al. Histomorphometric and fluorescence microscopic analysis of bone remodelling after installation of implants using an osteotome technique. *Clin Oral Implants Res*, 2002;13(6):595-602.
06. Frost HM. *Intermediary Organization of the Skeleton*. 1st ed. Boca Raton, FL: CRC Press; 1986:109-164.
07. Burri C, Wolter D. [The compressed autogenous spongiosis transplant (author's transl)]. *Unfallheilkunde*, 1977;80(5):169-175.
08. Hallidin A, et al. The effect of static bone strain on implant stability and bone remodeling. *Bone*, 2011;49(4):783-789.
09. Duncan RL, Turner CH. Mechanotransduction and the functional response of bone to mechanical strain. *Calcif Tissue Int*, 1995;57(5):344-358.
10. Kold S, et al. Compacted cancellous bone has a spring-back effect. *Acta Orthop Scand*, 2003;74(5):591-595.
11. Trisi P, et al. Implant micromotion is related to peak insertion torque and bone density. *Clin Oral Implants Res*, 2009;20(5):467-471.
12. Pagliani L, Sennerby L, Petersson A, et al. The relationship between resonance frequency analysis (RFA) and lateral displacement of dental implants: an in vitro study. *J Oral Rehabil*, 2013;40(3):221-227.
13. Trisi P, Colagiovanni M, Perfetti G. Implant Stability Quotient (ISQ) vs Direct In Vitro Measurement of Primary Stability (Micromotion): Effect of Bone Density and Insertion Torque. *Journal of Osteology and Biomaterials*, 2010;1(3).
14. Huwais S, Meyer EG. A Novel Osseous Densification Approach in Implant Osteotomy Preparation to Increase Biomechanical Primary Stability, Bone Mineral Density, and Bone-to-Implant Contact. *Int J Oral Maxillofac Implants* 2017;32:27-36
15. Lahens B, Neiva R, Tovar N, Alifrag AM, Jimbo R, Bonfante EA, et al. Biomechanical and histologic basis of osseodensification drilling for endosteal implant placement in low density bone. An experimental study in sheep. *J Mech Behav Biomed Mater* 2016;63:56-65.
16. Lopez CD, Alifrag AM, Torroni A, et al. Osseodensification for enhancement of spinal surgical hardware fixation. *Journal of the mechanical behavior of biomedical materials*, 2017;69:275-281
17. Trisi P, Berardini M, Falco A, Vulpiani MP. New osseodensification implant site preparation method to increase bone density in low-density bone: In vivo evaluation in sheep. *Implant Dent* 2016;25:24-31.
18. Huwais S. Fluted osteotome and surgical method for use. US2013/0004918. US Patent Application, 3 January 2013.
19. Huwais S. Autografting osteotome. WO2014/077920. World Intellectual Property Organization, 22 May 2014.

Connect With Us:



*Clinician experience and judgement should be used in conjunction with the Densifying Reference Guide recommendation and suggested use protocols.

*Clinician experience and judgement should be used in conjunction with the Densifying Reference Guide recommendation and suggested use protocols.



This Brand Knows Bone®

The Company Behind The Innovation

Versah® LLC was founded in April of 2014 by Dr. Salah Huwais, a practicing Periodontist in Jackson, Michigan. Dr. Huwais created Versah in order to provide dental implant clinicians with an improved and unique way to prepare osteotomies for implant placement. While placing implants over his 20 years of practice, Dr. Huwais consistently experienced the limitations of excavating bone using conventional drills when preparing osteotomies.


His goal was to create a new instrument and procedure to maintain healthy bone while preparing osteotomies rather than remove it to simply make space. Osseodensification is the procedure that led to the creation of Versah. The procedure is accomplished using the patented instruments Dr. Huwais invented and named the Densah® Burs.

Contrary to drilling away bone, the Densah Burs can be rotated in a non-cutting counterclockwise (CCW) direction at 800 – 1500 rpms and when coupled with irrigation, hydrodynamically densify bone through compaction autografting or Osseodensification.^{14, 18, 19} The result is a consistently cylindrical and densified osteotomy. Consistent osteotomies and densification are important to implant primary stability and to early loading.^{14, 15, 16, 17}




The History Of Osteotomy Preparation


Since the early days of dental implantology, osteotomies have been prepared using standard drills designed for use in industrial applications. These drill designs have proved to be functional for dental applications as implant success rates have been satisfactory over time but osteotomy preparation techniques have still been lacking for various reasons. Standard drill designs used in dental implantology are made to excavate bone to create room for the implant to be placed. The goal in implant placement is to achieve implant stability which is believed to increase implant osseointegration success. This is more important in recent days due to popular immediate loading protocols being implemented into treatment by many clinicians. Removing bone is contrary to achieving the primary stability desired. Standard drill designs, in twist or fluted shapes, cut bone effectively but typically do not produce a precise circumferential osteotomy. Osteotomies may become elongated and elliptical due to chatter of the drills. In these circumstances, the implant insertion torque is reduced leading to poor primary stability and potential lack of integration. Osteotomies drilled into narrow bone sites may produce dehiscence, buccally or lingually, which also reduces primary stability and will require an additional bone grafting procedure which adds cost and healing time to treatment.



1. Osteotomy created with Standard Drilling



2. Small buccal dehiscence after standard drilling




3. Large buccal dehiscence after standard drilling

Densah® Bur Technology


The Innovation That Makes Osseodensification Possible

Densah® Bur Technology: Instantly Familiar Yet Dramatically Different


The Densah® Bur uses advanced, patented technology designed to precisely cut bone in the clockwise (CW) direction and also to densify bone in a non-cutting counterclockwise (CCW) direction with copious irrigation. Operating at a counterclockwise speed it is able to preserve bone and gradually densify and expand the bone to prepare the osteotomy which enhances primary stability.¹⁴ The Densah® Bur features multiple flutes within a tapered geometry and is designed to produce a faster feed rate with minimum heat elevation.¹⁴ The chisel edge is narrow to minimize force thrust and employs a point angle designed to reduce chatter. Outward pressure coupled with irrigation at the point of contact creates a hydrodynamic compression wave so that bone is compacted laterally and apically by continuously rotating and concurrently forcibly advancing the Densah® Bur. The taper design allows the surgeon to instantly lift away from contact to enhance irrigation. This, along with real-time haptic feedback, makes the Densah® Bur intuitive for the surgeon to make an instant adjustment to the advancing force.^{18, 19}




Dual Mode Action
Rotating Counterclockwise the Flute Back Rake Angle Creates Osseodensification (Densifying Mode) Rotating Clockwise the Flute Edges Precisely Cut Bone (Cutting Mode)



Hydrodynamic Compression Wave
Facilitates Trans-Crestal Sinus Augmentation



Haptic Feedback
Real Time Haptic Feedback that is Intuitive to Inform the Surgeon if More or Less Force is Needed in an Instant^{18, 19}



Primary Stability
Compaction Autografting Facilitates Bone Springback Effect Which Enhances Implant Stability^{14, 15, 16, 17}

The Gold Tip™ Design
Tip Design, with the Flutes, Facilitates Compaction Autografting^{18, 19}

*Measure the drilling depth of the Densah® Bur from the widest part of its tip to the indication line. Regardless of the Densah® Bur diameter, the maximum additional tip depth is 1.0 mm. The accuracy of laser markings are tested within +/- .5mm.

Discover the Difference

Osteotomies Created With Standard Drills	Osteotomies Created With Densah® Bur Clockwise Cutting	Osteotomies Created With Densah Bur Counterclockwise Osseodensification
1. Osteotomy prepared in porcine Tibia with standard fluted drills	1. Osteotomy prepared in tibia with Densah Bur rotating clockwise (cutting)	1. Osteotomy prepared in tibia with Densah Burs rotating counterclockwise (Osseodensification)
2. uCT radiograph of the above osteotomy	2. uCT radiograph of the above osteotomy	2. uCT radiograph of the above osteotomy
3. Histology cross section of implant placed in osteotomy prepared using standard drills at Day 0	3. Histology cross section of implant placed in osteotomy prepared using Densah Burs rotating clockwise at Day 0	3. Histology cross section of implant placed in osteotomy prepared using Densah Burs rotating counterclockwise at Day 0
4. uCT radiograph of full length osteotomy prepared with standard drills	4. uCT radiograph of full length osteotomy prepared with Densah Burs rotating clockwise	4. uCT radiograph of full length osteotomy prepared with Densah Burs rotating counterclockwise



Scan this QR Code to view our Introduction Video

Spring Back Effect

When run in a non-cutting counterclockwise (CCW) direction with copious irrigation at 800-1500 RPM, a compression wave is produced, thereby auto-compacting the bone laterally and apically creating a spring back effect. This phenomena facilitates and enhances implant primary stability and increases BIC.¹⁴

Histological Cross-Section of Bone Healing with Osseodensification

Fig. 1-3: Six weeks healing in vivo, histometric results suggest that osseodensification positively influenced osseointegration.^{15, 16, 17}




Come Train with Us to Optimize Your Learning Curve

Osseodensification Hands-On Training

- Objectives:**
- How we can preserve bone, plastically expand it, and enhance its strength in a controlled method
 - The effect of Osseodensification on implant total stability development for early and immediate loading.
 - Risk assessment and avoiding complications.
 - Learn and practice the clinical versatility of Osseodensification.

Register Online at [versah.com/training/](https://www.versah.com/training/)

An Entire System of Minimally Invasive Procedures



"I attended your course in December in New York and am writing to tell you what a pleasure it is to be able to provide subantral augmentation and osseodensification this way. Thank you for your great innovation developing these instruments and technique. A real game changer. I have been doing implants for 25 years, and I think this is the best innovation in many years."
- Dr. Edwin Ginsberg, Periodontist. Great Neck, NY

"I attended the Milan course and since then have used your kit at least 4 times; always having the same great easy incredible result. You said that you'd change our clinical life... you did. How can I thank you? I'd like to let every dentist know how your burs are so incredible."
- Dr. Marco De Angelis, Oral Surgeon. Milan, Italy

"I studied the Versah videos while using Versah's drills for 6 months before taking the Hands-On course. The course was invaluable and cleared up a lot of fine-tuning questions I had about this amazing Osseodensification technique. Easily one of the best CE courses I've taken in a few years."
- Dr. Gregory Cumberland, General Dentist. Alberta, CAN

Hands-On Training Options

1-Day Hands-On Training Simulation model course 7 CE Credits	3-Day Live Surgery Simulation model & active participation in surgery 24 CE Credits*
---	---

The courses will be directed by **Dr. Salah Huwais**, the founder of Osseodensification, and Versah world-renowned faculty. Participants will practice on actual bone specimens and specially developed clinical simulation models. Dr. Salah Huwais maintains a private practice limited to Periodontics and Implantology. He is a Diplomat of The American Board of Periodontology and the American Board of Implantology; an adjunct assistant Professor at University of Minnesota. Dr. Huwais is the inventor of the Densah® Bur technology. He is the recipient of the 2017 AAID Foundation Isaih Lew Memorial Research Award.



Versah, LLC is designated as an Approved PACE Program Provider by the Academy of General Dentistry. The formal continuing education programs of this program provider are accepted by AGD for Fellowship, Mastership, and membership maintenance credit. Approval does not imply acceptance by a state or provincial board of dentistry or AGD endorsement. The current term of approval extends from 9/1/2017 to 8/31/2020. Provider ID# 371646

*7 CE credits from Versah for Hands-On Simulation Course.
17 CE credits from Augusta University Dental School for LIVE-Surgery Course.