

메가젠이 추천하는

세계 초일류 상품

WORLD NO. 1

Anatomic Healing Abutment™

치아 위치별 이상적이고 심미적인 Emergence Profile을 형성을 위한 Solution




# Anatomic Healing Abutment™

- 자연치와 같은 심미적인 최종보철 제작
- 발치와에 딱 맞는 형상
- 즉시사용 가능! 체어타임 감소

## 최적의 잇몸 Healing! Save time, Save money


- 치아 형상에 맞는 디자인으로 위치 마다 형상을 다양화하여 발치와 형상에 맞는 제품을 선택가능  
잇몸을 자연치가 있을때와 같이 유지
- 레진을 이용하여 직접 제작하거나 밀링 머신을 사용하지 않고 기성품 선택할 수 있어 제작에 소요되는 시간과 비용을 save!



[일반 Healing 사용시]



[Anatomic Healing 사용시]



## 자연치와 같은 심미적인 최종보철 제작

- 자연치 그대로의 잇몸을 형성하여 심미적인 최종보철 제작할 수 있도록 도와주는 역할
- 임상 의와 환자 모두가 만족스러운 결과를 얻을 수 있음



Courtesy of Dr. Ifan Abas

## Digital Abutment와 최적의 콜라보레이션

- Anatomic Healing Abutment로 자연치가 있을때의 잇몸을 유지시키고 여기에 Customized Abutment로 최종보철을 제작하여 환자 맞춤형 보철을 완성

### Clinical Cases

Courtesy of Dr. Sam Omar

Case 1.



Case 2.



# Anatomic Healing Abutment™

All in one package



Order code

System	Height(mm)	Type	Ref.c
ANYRIDGE® by MEGAGEN	4	Hex	C-AHAR04HP
		Non-Hex	C-AHAR04NP
	5	Hex	C-AHAR05HP
		Non-Hex	C-AHAR05NP
	7	Hex	C-AHAR07HP
		Non-Hex	C-AHAR07NP
BLUEDIAMOND IMPLANT	4	Non-Octa	C-AHBD04NP
		Octa	C-AHBD04OP
	5	Non-Octa	C-AHBD05NP
		Octa	C-AHBD05OP
	7	Non-Octa	C-AHBD07NP
		Octa	C-AHBD07OP
ANYONE® by MEGAGEN	4	Hex	C-AHAO04HP
		Non-Hex	C-AHAO04NP
	5	Hex	C-AHAO05HP
		Non-Hex	C-AHAO05NP
	7	Hex	C-AHAO07HP
		Non-Hex	C-AHAO07NP

## Incisor



Incisor				
MD (mm)	4.0	4.5	6.0	7.0
LL (mm)	5.0	4.5	5.0	6.0
Position No.	I1	I2	I3	I4

## Canine



Canine	
MD (mm)	5.0
LB (mm)	5.5
Position No.	C1

## Pre-Molar



Pre-Molar		
MD (mm)	4.5	5.0
LB (mm)	6.0	7.0
Position No.	P1	P2

## Molar



Molar								
MD (mm)	6.0	6.0	6.0	7.0	7.0	7.0	8.0	8.0
LB (mm)	7.0	8.0	9.0	8.0	9.0	10.0	9.0	10.0
Position No.	M1	M2	M3	M4	M5	M6	M7	M8

## Special



Special											
MD (mm)	4.5	5.0	5.0	6.0	6.0	6.0	7.0	7.0	7.0	8.0	8.0
LB (mm)	6.0	6.5	7.0	7.0	8.0	9.0	8.0	9.0	10.0	9.0	10.0
Position No.	S1	S2	S3	S4	S5	S6	S7	S8	S9	S10	S11

REV.03