

Titanium Brush



Best item for Peri-implantitis

Peri-implantitis

I. Titanium Brush / Set

- Excellent for removing bacteria and other substances adhered to exposed implants.
- Highly effective for treating peri-implantitis and peri-coronitis without removing the prosthetic.

Description	Ref.C
Titanium Brush (Long)	DTST-6352
Titanium Brush (Short)	DTST-6353
Titanium Brush Set (10pcs)	DTST-TB-SET

1. Features

- Made of special materials to minimize implant damage
- Optimized design for maximum accessibility

2. Set Component




Set: Long(5ea), Short(5ea)

3. Single Items



4. Test Case



5. How to Use



Recommended RPM: 800~1200

Titanium Brush



Head Office & Factory
Gangnam Office

45, Secheon-ro 7-gil, Dasa-eup, Dalseong-gun, Daegu, Korea
MegaGen Tower, 607 Seolleung-ro, Gangnam-gu, Seoul, Korea

T. +82-1544-2285
T. +82-1566-2338