

EASY EXTRACTION KIT

by MEGA'GEN





Easy Extraction KIT

005 I. Easy Extraction KIT

006 1. EASY PERIOTOME

008 2. EASY LUXAOR

010 3. EASY FORCEP

011 4. SURGICAL CURETTES



AnyRidge is Perfect for Immediate implant placement following extraction with Easy Extraction Kit

Easy & atraumatic extraction, possible with Easy extraction Kit.

EASY PERIOTOME


“AMAZING FLEXIBLE”



EASY PERIOTOME adopts the flexible type which can be bent freely, and ensures approach to the tooth root, so that you can cut off the periodontal ligament very easily.

EASY LUXATOR


“ELASTIC & THIN”



EASY LUXATOR is made of a special alloy, so it has elasticity and is thinner than other Luxators. It has the advantage of free from breakage, preserving the shape of the skeleton, and minimizing tissue damage.


EASY FORCEP

“DECREASE the possibility of ROOT FRACTURE”



EASY FORCEP delivers force by grasping more apical portion than previous instruments, decreasing the possibility of root fracture significantly during extraction.

Non-invasive exodontia & Biology-Driven Implant System will make better esthetics!



- Courtesy of Dr. Achraf Souayah

Have you made the PARADIGM SHIFT yet?

Do it the AnyRidge Way

Excellent primary stability with Knifethread®

- For smooth insertion & stronger primary stability
- No cutting edge for minimal invasion
- Perfect results with any bone condition

Excellent secondary stability with Xpeed®

- Faster & stronger osseointegration
- Surface technology incorporating Ca²⁺ ions on S-L-A treated surface
- 100% acid-residue-free surface



I. Innovation of Exodontia **Easy Extraction KIT**

Non-invasive exodontia will make better esthetics!

No more difficult & complicated extraction!!

Easy & atraumatic extraction, possible with Easy extraction Kit.

EASY Extraction KIT Basic set

(3306, 3308, 3251, 3252, 3113, 3114)

Ref.C

3000



EASY EXTRACTION BASIC KIT IS USEFUL FOR THE FOLLOWING CASES

- **Immediate implant placement following extraction**
- Extraction of root canal treated teeth
- Intentional tooth replantation
- Residual root with high possibility of fracture

Upgrading conventional areas through matchless design

Unlike other existing extraction tools, 'Easy Extraction Kit' has been developed in order to provide more safe and precise extraction surgery while considering the convenience of the user.

MegaGen highly recommends using 'Easy Extraction Kit' to minimize the damage of the alveolar bone during tooth extraction and calcification.



I. Easy Extraction KIT. Easy Extraction KIT

1. EASY PERIOTOME

This product is used for minimizing damage to the alveolar bone as it helps cut the periodontal ligament to detach teeth from the periodontal ligament to place implant immediately after tooth extraction or general tooth extraction.

Characteristics

FLEXIBLE TYPE

EASY PERIOTOME adopts the flexible type which can be bent freely, and ensures approach to the tooth root, so that you can cut off the periodontal ligament very easily.



Various Products

2mm STRAIGHT




3mm STRAIGHT



2mm CURVED



3mm CURVED



2mm CURVED



2mm CURVED





I. Easy Extraction KIT

2. EASY LUXATOR

It is a New Easy Luxator which has flexible and thin blade unlike other Luxators. Due to this structural feature, the tip of the Luxator easily digs into the PDL space. So the distortion of teeth from the alveolar bone can be induced with even a small force. It also enables Atraumatic extraction because it minimizes the damage of surrounding soft tissue and alveolar bone.

Characteristics


It is made of a special alloy, so it has elasticity and is thinner than other Luxators. It has the advantage of free from breakage, preserving the shape of the skeleton, and minimizing tissue damage.



EASY LUXATOR is used for cutting off only the periodontal ligament to make easier extraction. Do not lift it up or twist with excessive force such as ELEVATOR

Various Products


3mm STRAIGHT



Ref.C

3111
3112

4mm STRAIGHT



Ref.C
3113
3114

4mm CURVED



Ref.C
3115
3116


I. Easy Extraction KIT

3. EASY FORCEP

'EASY FORCEP' provides a new specialized tip design, not the existing way to hold a tooth for extraction, but a new approach to penetrate into the inside of periodontal tissue for tooth extraction. It lessens the luxation effect and possibility of root residues while minimizing damage to the alveolar bone.

Characteristics

EASY FORCEP delivers force by grasping more apical portion than previous instruments, decreasing the possibility of root fracture significantly during extraction.




Various Products

Maxillary Anterior tooth



Ref.C
3251



Mandibular Anterior tooth



Ref.C
3252





Mandibular Posterior tooth



Ref.C

3253



Root [UNIVERSAL]



Ref.C

3254




IMPLANT



Ref.C

3220



I. Easy Extraction KIT

4. Surgical Curettes

Product




Ref.C

4004

Fine Saw Tip

Saw type tip allows easy curettage of soft tissue.



EASY EXTRACTION KIT

by MEGA[!]GEN

